




Data Sheet

S-60-I Medical Device



The **S-60-I** medical device, composed of an electronic board (prefix **MB**), a micromotor (prefix **MO**) and a cable connection (prefix **HO**), is intended to drive a straight or contra-angle handpiece as part of a dental treatment unit.

The **S-60-I** medical device is intended for the following medical indications:

- Implantology: drilling holes into bone tissue
- Prophylaxis: hygienist, polishing, etc.

This equipment is intended for professional use and requires installation from qualified personal.

MO-94 motor

The motor converts electrical energy to rotation and transfers rotation to the instrument tool. This brushless motor is able to maintain any speed between 0 and 40 000 rpm with controlled torque and is equipped with optional air cooling and powerful white lighting. Dassym motors withstand thousands cycles of sterilization.

This motor is especially designed to drive an ISO 14457 & ISO 3964 compliant accessory such as a straight or contra-angle handpiece to which a rotating or oscillating tool is attached. Conjugated to the precision of the MB-60 closed loop three phases sine drive, the MO-94 motor is particularly suited for high torque at low speed operation.

The external parts of the motors are made of stainless steel (X10CrNiS18).

Technical specifications

Dimensions	MO-94-LA	Unit
Motor length excluding coupling	61	mm
Motor maximum diameter	21.8	mm
Performances	MO-94-LA	Unit
Speed range	0 – 40 000	rpm
Maximum torque	5	Ncm
Maximum power	80	W _{MEC}
Characteristics	MO-94-LA	Unit
Optional cooling air pressure	300	kPa
Cooling air flow consumption	< 40	l/min
Cooling air flow at the nozzle	—	
Features	MO-94-LA	
Sterilizable without protection	✓	
Full torque available from standstill	✓	
Rocking & auto reverse modes	✓	
Powerful white LED lighting	✓	
Optional motor air cooling	✓	
ISO 3964 & INTRAmatic Lux® coupling	✓	

INTRAmatic Lux® is a Kaltenbach & Voigt, Biberacht, DE, trade mark.

HO-05 cable connection

The cable conveys electrical energy from the dental unit to the motor. It allows disconnection of the motor side with an optional unit side detachable connector for sterilization.

The connection cable comprises a Dassym proprietary motor connector with its nut and encloses a length of five wires cabling for dental unit connection.

The external parts of the cable are made of silicone and stainless steel (X10CrNiS18).

Technical specifications

Dimensions	HO-35	Unit
Cable length	1.6	m
Cable mean diameter	5.3	mm
Features	HO-35	
Stranded wires for motor phases and lighting	✓	
Optional unit side connector for sterilization	✓	

MB-60 driver board

The motor driver board receives commands from the dental unit and sends controlled electrical energy to the motor. The board is to be integrated in the dental unit.

The board monitors its working conditions and finely drives the motor speed, torque and lighting. The motor is safeguarded by electronic limitation of its temperature.

Technical specifications

Dimensions	MB-60			Unit
Length × Width × Height	102 × 59 × 25			mm
Fastening holes distance & diameter	95 × 51 Ø 3.2			mm
Characteristics	Min	Nominal	Max	Unit
Supply Voltage	24	32 – 36	48	V _{DC}
Maximum Consumption	180	240	360	VA
Features	MB-60			
Digital protocols & analog controls	✓			
Optional pressure speed reference control	✓			
Optional router for two motors connection	✓			

Main features

- brushless motors of the sturdiest construction manufactured in the Swiss watchmaking area
- speed range from standstill to 40 000 rpm with finely controlled torque up to 5 Ncm
- full torque from standstill with patented technology, rocking & auto reverse modes
- powerful LED white lighting compatible with INTRAmatic Lux® handpieces
- light, flexible and strong cable with proprietary connector for reliable motor connection
- highly configurable electronic controller board with versatile analog and digital interfaces
- ready to be installed as a replacement for an air turbine or a brushed DC electric motor
- direct and compatible replacement for previous Dassym MB-30, MB-92 & MB-92-B products
- extensive compliance with the highest international safety and performance standards
- classe IIa medical device in full conformity to the Regulation (EU) 2017 / 745 (MDR)

Orderable parts

UDI-DI	Trade name	Trade mark	Product reference	Composition	Part reference
EDAS1506004010	S-60-I	Dassym	150.60.941	MO-94-LA	102.08.200
				MB-60	251.60.000
				HO-05	121.05.001

For other options, please contact Dassym.



Dassym SA
Route de Séprais 19a
2856 Boécourt
SWITZERLAND

SRN: CH-MF-000048199
www.dassym.com
info@dassym.com

Active patents
CH709334, EP2858589



Veranex Germany GmbH
Landsberger Straße 302
80687 München
GERMANY

SRN: DE-AR-000005578
www.veranex.com
info@veranex.com

